

### General

The three axis HUBER in-vacuum polarisation analyser has been designed to facilitate the study of changes in the polarisation of the x-ray beam after diffraction from the sample.

The in-vacuum design allows the device to be used at low energies (ca. 3keV).

### Technical details

The polarisation analyzer comes with a T-shaped vacuum chamber, clamped on a CF vacuum flange with a support for the crystal holder. The top flange with a pumping connection can be removed in order to get access to the analyser crystal.

On the  $2\theta$  linear translation, which seals the chamber with a Kapton window, a detector with a max. diameter of 48mm with shielding can be mounted. The linear range is  $\pm 20$ mm which corresponds to  $2\theta = \pm 7^\circ$ .

Optionally a further detector can be mounted in direction of the scattered beam on the back side of the support of the of the  $\sigma/\pi$ -axis (Huber model 410a) and can be used windowless or with Kapton window in vacuum for measurements without analyzer crystal.

The whole assembly is mounted on a prism slide (or optional other slides) to be fitted on the  $2\theta$  arm of a diffractometer.

### Specification

Ranges [°]:

$\theta$  (model 409) [°]:  $\pm 5$

$2\theta$  (linear translation at  $90^\circ$  fixed angle [°]:  $\pm 7$

$\pi/\sigma$  (model 410a) [°]:  $\pm 180$


ca.  $\pm 110$   
(with detector mounted)

Resolution [°]: 0,005 to 0,0001  
(depending on application and motorisation)

Vacuum chamber with 2 ports for the direct beam with DN40KF flanges

One port for vacuum pumping

Mounting flange for analyzer Crystal



### Applications

The  $\Theta$ -axis (Huber model 409) allows alignment of the analyser crystal to a diffracting condition, with  $\Theta$  Bragg at, or close to,  $45^\circ$ . The second rotation axis (Huber model 410a) allows rotation of the diffracting plane of the analyser crystal about the beam.

Any component of the incident polarisation lying in the diffracting plane will go to zero on charge scattering for  $\Theta$  Bragg =  $45^\circ$ . The third linear axis allows translation of the detector in a  $2\Theta$  geometry in order to track the diffracted beam.

Conventionally, for vertical scattering experiments, the incident polarisation is referred to as  $\sigma$  polarised and any component orthogonal to it (i.e. vertically polarised) as  $\pi$  polarised. Thus, by positioning the rotation about the beam such that the diffracting plane of the analysing crystal is vertical, one is sensitive only to the  $\sigma$  polarised component and conversely a horizontal analyser crystal diffraction plane measures the  $\pi$  component.

Any additional components of polarisation introduced on diffraction from the sample, due to its magnetism, may therefore be quantified.

