

Reference	3496	Description of Space or Activity/Task or Equipment	Use of the Horiba LabRam HR Spectrometer for Raman and photoluminescence measurements conducted at room temperature, located in Millburn house G78
Assessment Date	06/09/2024	Publish To Portal	No
Assessor Name	Ben Breeze	Risk Assessment Title	Horiba LabRam evo HR
Assessment Team Members		Review Date	No Review Set
Role / Space / Project Reference		Current Risk Level (1=Very Low, 2=Low, 3=Moderate, 4=High, 5=Very High)	2
Department	Use the search function above or double click here for org chart -> Academic Faculties -> Faculty of Science, Engineering and Medicine -> Research Technology Platforms - RTPs -> Spectroscopy Research Technology Platform - RTP	Final Risk Level (1=Very Low, 2=Low, 3=Moderate, 4=High, 5=Very High)	2
Location Details	Central Campus-Millburn House-Ground Floor-G78 - (01.005.000.051) Spectroscopy RTP Laboratory Millburn House G78	Risk Assessment Number	0



Risk Assessment Category				Additional Information			<ul> <li>Users are responsible for having performed a risk assessment for their samples which must be logger at:</li> <li>https://warwick.ac.uk/research/rtp/spectroscopy/Sa ety/Sample</li> <li>When booking all users will provide a reference to the risk assessment for the sample used.</li> <li>All Users must complete training course before operating the spectrometer.</li> <li>Bookings are restricted to trained users who are recorded in the px-pl-users web group.</li> <li>All users will require retraining if they have not operated the system for 18 months.</li> <li>Document History</li> <li>Version Date Reviewer Comments</li> <li>1 05/01/18 Ben Breeze Re reviewed when Equipment was brought under RTP Control</li> <li>2 24/01/19 Ben Breeze</li> <li>3 08/09/21 Ben Breeze</li> <li>Scheduled Review Updated Format updated, minor changes to user groups and description</li> </ul>				
Date Record Create	d 06/09/2024										
Hazard Type & Hazard Description	Who may be at Risk? & How May Person(s) Be Harmed	Existing Control Measures	L	S	R	VH, what ad	nt risk is M, H or Iditional Control are required?	L	S	R	
NIR Class 3B/4 laser or Class 1 laser product Laser strike (burns, loss of sight) 488 nm 660 nm	Laboratory Worker Staff Student Light affecting eye or skin	Admin - Refer to more detailed risk assessment for Class 3B/4 lasers or for Class 1 laser products.	Serious	Unlikely	Low			Serious	Unlikely	Low	



Equipment/Plant Trapping or object being dropped. Contact with rotating mechanical parts (Friction). Contact with mechanical pinch or nip points.	Laboratory Worker Staff Student Trapped fingers Bruies	Admin - Information, instruction, supervision & training. Eng/Admin - Adjustable guarding in place to allow safe use. Admin - Awareness training provided. Admin - Controls are suitably located and identified. Admin - Safe System of Work (SSoW) established and trained out.	Minor	Possible	Low		Minor	Possible	Low
Equipment/Plant broken glassware Contact with sharp blade (cuts).	Laboratory Worker Staff Student Cuts from broken quartz	Admin - Defective equipment taken out of use. Admin - Information, instruction, supervision & training. PPE - Personal Protective Equipment (PPE) is issued and worn (gloves goggles). Take care, dispose of sharps and broken glass in the approved receptacles provided.	Minor	Unlikely	Very Low		Minor	Unlikely	Very Low
Electricity Contact with live electrics.	Laboratory Worker Staff Student Electric Shock and/or burns	Eng - Additional power sockets installed (reducing need to daisy-chain extension cables). Eng - Cables and leads are appropriately insulated. Eng - Fixed guarding is in place preventing access.	Serious	Unlikely	Low		Serious	Unlikely	Low
Assessment Conclusion Safe operating procedure									