

## Gloves Selection Guide




Glove Material	Uses
<b>Nitrile</b> <i>Disposable and re-useable</i>	<b>Good</b> with oils, greases, some weak acids and bases. <b>Avoid</b> intentional contact with ketones, oxidising acids and organic compounds containing nitrogen Good all round chemical splash protection
<b>Vinyl (PVC)</b>	<b>Good</b> for weak acids and bases, oils, fats, amines and peroxides <b>Avoid</b> intentional contact with ketones and aromatic solvents
<b>Butyl Rubber</b> <i>Re-useable</i>	<b>Good</b> for ketones and esters. <b>Poor</b> for petrol, aliphatic, aromatic and halogenated solvents and where dexterity is critical
<b>Polyvinyl alcohol (PVA)</b> <i>Re-useable</i>	<b>Good</b> for aromatics and chlorinated solvents <b>Avoid</b> for water based solutions (dissolves in water)
<b>Neoprene</b> <i>Disposable and re-useable</i>	<b>Good</b> for weak acid, strong bases, alcohols, fuels, peroxides, hydrocarbons and phenols. Remains flexible at low temperatures. <b>Poor</b> for halogenated and aromatic solvents
<b>Viton</b> <i>Re-useable</i>	<b>Good</b> for chlorinated and aromatic solvents <b>Poor</b> for ketones and dexterity
<b>Latex (Natural rubber)*</b> <i>Disposable and re-useable</i>	<b>Good</b> for incidental contact with phenol, dilute aqueous acids and bases, inorganic chemicals <b>Poor</b> for oils, greases and hydrocarbon derivatives

\*Natural rubber latex can cause a very serious allergic reaction. Always look for an alternative where possible. **Powdered latex gloves must not be used. If unpowdered gloves are used, the maximum free protein content should be 50 ug per gram of glove material.**

## Guidance for the use of gloves

- Select the gloves that is most appropriate for the chemicals you are handling
  - **Disposable** gloves are for small incidental splash or contact
  - **Re-useable** gloves are better suited when immersion is unavoidable
- Check that the glove is fit for purpose – no tears or holes
- Ensure the glove fits correctly
- Do not re-use disposable gloves
- Keep surface of glove clean when working – it may be necessary to wash the external surface of the glove regularly
- Remove gloves in the correct manner – see ‘Skin Care Guide’
- Wash hands after removing gloves
- Always remove gloves before :-
  - leaving the laboratory;
  - using taps, phones and keyboards;
  - writing in lab book
- Dispose of gloves in a proper manner
- If you contaminate your disposable glove with toxic chemicals change immediately



Check the packaging of gloves for **BS EN 374** and the symbols above and also PPE cat I – for minimal risk, PPE cat III for serious harm (CE marked)

Examination gloves are not PPE – they are medical devices.

Gloves tested to ASTM F1671 are resistant to viruses – particles >27nm.

**There is no single glove material that protects against all chemical hazards**

**Contact your Health and Safety Adviser for further guidance on the correct selection of gloves**

