

# Non-Traditional Biomaterial Solutions to Infection. Detect, Neutralize and Kill

The Gibson Group: Chemistry and Medical School

WARWICK  
INTEGRATE  
ANTIMICROBIAL RESISTANCE


M.I.Gibson@warwick.ac.uk

@LabGibson




www2.warwick.ac.uk/go/gibsongroup

## 1. Pathogen Detection

### Gold nanoparticles as colorimetric sensors




- Gold nanoparticles (AuNPs) have interesting optical properties. Red in solution → blue upon aggregation.
- Glycosylated AuNPs (glycoAuNPs) could be used to detect bacteria and lectins.



*J. Mater. Chem. B.*, 2014, 2, 1490; *ACS Macro Lett.*, 2014, 3, 1004;  
*J. Mater. Chem. B.*, 2016, 4, 3046; *Analyst*, 2017, 141, 4305


### High density glycan arrays

- Glycan terminated polymers immobilised onto surface functionalised glass
- Via Direct microcontact printing
- Generation of high-density arrays of surface immobilised glycans, with polymer tethers for reduced non-specific binding
- Multiplexing of data enables 'barcoding' of pathogen binding patterns to enable discrimination




*Biomater. Sci.*, 2015, 3 (1), 175-181; *RSC Adv.*, 2015, 5, 53911-53914;  
*Mol. BioSyst.*, 2016, 12, 341-344; *Biomacromolecules*, 2016, 17(8), 2626-2633

## 2. Neutralize





- Anti-adhesion therapy prevents binding of bacteria and toxins to host cells
- Does not target bacteria directly, so less evolutionary stress
- Reduced the need for antibiotics
- Glycopolymers have high affinity due to presentation of 10's of copies of each sugar
- Here we engineer high-affinity binders without total synthesis
- Demonstrated prevent binding of cholera-toxin to GM-1, stopping infection



*Angew. Chem. Int. Ed.*, 2012, 51, 7812;  
*Chem. Sci.*, 2014, 5, 1611

## 3. Kill

- Using cationic polymers as non-traditional antibiotics
- Not previously evaluated against mycobacteria, which includes Tuberculosis
- We found that mycobacteria were killed selectively over gram negative bacteria
- Nature of cationic polymer was crucial for selectivity



- Mechanism of action
- Polymers were found not to lyse the bacterial membrane
  - Electron microscopy suggested stress of the cell walls
  - Confocal microscopy did not indicate membrane permeation

*Biomacromolecules*, 2017, 18, 1592

Collaborators  
(Warwick)  
Dr E. Fullam  
Dr N. Waterfield

Iceni  
Diagnostics

The Leverhulme Trust

BBSRC  
bioscience for the future

erc  
European Research Council

EPSRC  
Engineering and Physical Sciences  
Research Council

MRC  
Medical Research  
Council